


VacuDAP goes wireless

Install the Dose Area Product meter VacuDAP at your X-ray machine and save the interface cable by using  **Bluetooth**[®] wireless technology



- * save expenses for cable installation between DAP meter and operating desk
- * assure complete functionality of the standard VacuDAP
- * display DAP and DAP rate simultaneously

Specifications:

Bluetooth type	Class 1/ +7 dBm (5 mW)
Antenna	internal
Operating range	> 20 m
Operating mode	point-to-point communication
Output frequency	2.402 - 2.480 GHZ, ISM band radio type approved for Europe, US, CA

VacuTec Meßtechnik GmbH
Dornblüthstraße 14
01277 Dresden

Tel.: +49 (0)351 31724-0
Fax.: +49 (0)351 3105085
e-mail: info@vacutec-gmbh.de
www.vacutec-gmbh.de